



Cardinal Health™

Protexis™ Surgical Gloves





Contents

3 Introduction

4-8 Protexis™ Surgical Gloves Perfect Fit

9-23 Protexis™ Surgical Gloves Portfolio

NON-LATEX POLYISOPRENE

- 10 Protexis™ PI
- 11 Protexis™ PI with Neu-Thera™
- 12 Protexis™ PI Blue with Neu-Thera™
- 13 Protexis™ PI Micro
- 14 Protexis™ PI Orthopaedic
- 15 Protexis™ PI Ortho

NON-LATEX NEOPRENE

- 16 Protexis™ Neoprene
- 17 Protexis™ Neoprene Essential

LATEX

- 18 Protexis™ Latex
- 19 Protexis™ Latex with Neu-Thera™
- 20 Protexis™ Latex Essential
- 21 Protexis™ Latex Blue with Neu-Thera™
- 22 Protexis™ Latex Micro
- 23 Protexis™ Latex Ortho

24 Environmental Stewardship

25-29 Testing Standards and Data

30-34 Product Quick Reference

Cardinal Health™ Protexis™ Surgical Gloves

Protecting the hands that heal.



You give everything to your patients, it's who you are. Just as it's in your instinct to provide care, it's ours to contribute to your safety every step of the way. To protect you and your hands that heal.

That's why we build quality, clinical best practice, and trust directly into every Protexis™ Surgical Glove, to provide you with the **Perfect Fit** for every type of procedure.



Quality




Clinical best practice



Trust

Are you ready to find your Perfect Fit?



QUALITY

Protexis™ Surgical Gloves

Perfect Fit



CLINICAL BEST PRACTICE



TRUST

Protexis™ Surgical Gloves

Perfect Fit



QUALITY



CLINICAL BEST PRACTICE



TRUST

Quality

You shouldn't have to sacrifice quality for a comfortable fit and feel. That's why for more than five decades we've been committed to crafting proprietary molds and formulations that mimic the actual anatomy and physiology of your hands and owning the entire manufacturing and supply chain process to ensure the highest quality standards are adhered to at every step.



2,000 quality variables monitored



AQL of **<0.65**
exceeding industry standards



ISO, EN, ASTM and PPE
standards met globally



"The gloves fit well and
are comfortable."

RN, Top 3 *United States News & World
Report* Hospital¹



Protexis™ Surgical Gloves

Perfect Fit



QUALITY



CLINICAL BEST PRACTICE



TRUST

A mold unlike any other.

There are important differences in surgical glove manufacturing and technology that can impact the way a glove feels and performs. We believe that the surgical glove mold used makes all the difference in the fit and feel of the glove.

The Cardinal Health™ Surgical Glove mold features an independent thumb design — meaning the finger placement replicates the anatomical position of a resting hand — allowing for natural, comfortable movement.


One study shows that gloves featuring an independent thumb require less force for similar displacements of the thumb than those gloves produced without an independent thumb feature.¹

No matter what type of surgery you are in, surgical gloves produced on a mold with an independent thumb design ultimately can help reduce hand fatigue so you can perform at your best.

Less flexion and extension force required across the palm due to independent thumb design.

Anatomic thumb to reduce thumb and hand fatigue.

Interlocking beaded cuff to help prevent roll down.



1. Bennet, M. A. & Tekamp, D. A. Surgical Glove Comparison Testing. Stress Engineering Services. February 2012.

Protexis™ Surgical Gloves Perfect Fit



QUALITY



CLINICAL BEST PRACTICE



TRUST

Clinical best practice

The Perfect Fit is a commitment to provide your team with the most clinically recommended product options that meet their unique needs, as well as the ongoing education and tools to maximize their safety in the operating room.



“Education and consistency during our conversion of all gloves with Cardinal Health helped to garner staff and physician buy-in.”

Executive Director, 70+ facility health system in the United States¹



Videos



Clinically-authored white papers and journals



Posters, fact sheets and interactive tools



Sampling, sizing and product evaluation support

TOPICS INCLUDE:

Latex safety

Double-gloving

Dermatitis prevention



Protexis™ Surgical Gloves

Perfect Fit



QUALITY



CLINICAL BEST PRACTICE



TRUST



Trust

The Perfect Fit is a commitment to providing support throughout the conversion process to reduce the burden of change for your entire team. We have developed best-in-class conversion processes to help you improve efficiency, while reducing inventory and costs through standardization.

In fact, **79%** of our customers in the United States **actually saved money** when they went latex-free with Protexis™ Surgical Gloves.²



"The Cardinal Health team made the conversion seamless."

VP of Perioperative Services,
Large Academic Teaching Institution
in the United States¹

1. GHX Units, All Channels, 2017.

2. Techvalidate Research Study, December 2017.



Cardinal Health™ Protexis™ Surgical Gloves Portfolio

Protexis™ Surgical Gloves

NON-LATEX POLYISOPRENE

- 10 Protexis™ PI
- 11 Protexis™ PI with Neu-Thera™
- 12 Protexis™ PI Blue with Neu-Thera™
- 13 Protexis™ PI Micro
- 14 Protexis™ PI Orthopaedic
- 15 Protexis™ PI Ortho

NON-LATEX NEOPRENE

- 16 Protexis™ Neoprene
- 17 Protexis™ Neoprene Essential

LATEX


- 18 Protexis™ Latex
- 19 Protexis™ Latex with Neu-Thera™
- 20 Protexis™ Latex Essential
- 21 Protexis™ Latex Blue with Neu-Thera™
- 22 Protexis™ Latex Micro
- 23 Protexis™ Latex Ortho



This product is not made
with natural rubber latex

Protexis™ PI Surgical Gloves

- Engineered to protect in a wide array of clinical cases
- A multipurpose solution that offers tactile response with barrier protection
- Synthetic, not made with natural rubber latex

Catalog no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D72PT55X	5.5	287mm	0.23mm	0.17mm	0.17mm	Synthetic polyisoprene (PI)	 Cream	Beaded/ Rolled	50	200
2D72PT60X	6									
2D72PT65X	6.5									
2D72PT70X	7	300mm								
2D72PT75X	7.5									
2D72PT80X	8									
2D72PT85X	8.5									
2D72PT90X	9									



See Appendix page 29 for complete testing standards



See Appendix page 28 for chemotherapy agent permeation results






This product is not made with natural rubber latex

Protexis™ PI with Neu-Thera™ Surgical Gloves

- Designed to be comfortable and reliable for a broad range of surgical procedures
- Same great engineering as our Protexis™ PI Surgical gloves with our patented **Neu-Thera™** Coating
- Synthetic, not made with natural rubber latex

⦿ What is Neu-Thera™ Coating?¹

We place moisturizing Neu-Thera™ Coating on the inside of Protexis™ PI with Neu-Thera™ Surgical Gloves to **help promote overall skin wellness and support hand health.** It provides relief of dry, flaky skin and minimizes skin shedding.

Catalog no.	Size	Length	Thickness ²			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D73TE55	5.5	287mm	0.23mm	0.25mm	0.17mm	Synthetic polyisoprene (PI) with Neu-Thera™ Emollient Coating	 Cream	Beaded/ Rolled	50	200
2D73TE60	6									
2D73TE65	6.5									
2D73TE70	7	300mm								
2D73TE75	7.5									
2D73TE80	8									
2D73TE85	8.5									
2D73TE90	9									



See Appendix page 29 for complete testing standards

1. Data on file
2. Thickness tested in accordance with ASTM D 3577




This product is not made
with natural rubber latex

Protexis™ PI Blue with Neu-Thera™ Surgical Gloves

- Good underglove solution when double-gloving
- Distinct blue color aids in alerting wearers to perforations in the outer glove
- Synthetic, not made with natural rubber latex

⑤ What is Neu-Thera™ Coating?¹

We place moisturizing Neu-Thera™ Coating on the inside of Protexis™ PI Blue with Neu-Thera™ Surgical Gloves to **help promote overall skin wellness and support hand health**. It provides relief of dry, flaky skin and minimizes skin shedding.

Catalog no.	Size	Length	Thickness ²			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D73EB55	5.5	287mm	0.20mm	0.14mm	0.14mm	Synthetic polyisoprene (PI) with Neu-Thera™ Emollient Coating	 Blue	Beaded/ Rolled	50	200
2D73EB60	6									
2D73EB65	6.5									
2D73EB70	7	300mm								
2D73EB75	7.5									
2D73EB80	8									
2D73EB85	8.5									
2D73EB90	9									



See Appendix page 29 for complete testing standards






This product is not made
with natural rubber latex

Protexis™ PI Micro Surgical Gloves

- Stretches and conforms to your hand contour, staying comfortably in place
- Good solution in a thin double-gloving system where fingertip sensation is essential
- Heightened tactile response with a comfortable, smooth, anti-slip finish
- Synthetic, not made with natural rubber latex

Catalog no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D73PM55	5.5	287mm	0.20mm	0.14mm	0.14mm	Synthetic polyisoprene (PI)	 Cream	Beaded/ Rolled	50	200
2D73PM60	6									
2D73PM65	6.5									
2D73PM70	7	300mm								
2D73PM75	7.5									
2D73PM80	8									
2D73PM85	8.5									
2D73PM90	9									



See Appendix page 29 for complete testing standards



See Appendix page 28 for chemotherapy agent permeation results






This product is not made with natural rubber latex



Protexis™ PI Orthopaedic Surgical Gloves

- Thickest glove in the synthetic portfolio
- Smooth finish for tactile sensitivity
- Water-based hydrogel coating for easy donnability
- Durable for broaching and tactile for pinning
- Rich brown color reduces glare from surgical lighting
- Synthetic, not made with natural rubber latex

Catalog no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D73HT60	6	291mm	0.34mm	0.26mm	0.21mm	Synthetic polyisoprene (PI) with water-based hydrogel polymer coating	 Brown	Beaded/ Rolled	40	160
2D73HT65	6.5									
2D73HT70	7									
2D73HT75	7.5									
2D73HT80	8									
2D73HT85	8.5									
2D73HT90	9	305mm								



See Appendix page 29 for complete testing standards






This product is not made with natural rubber latex

Protexis™ PI Ortho Surgical Gloves

- Designed for use in orthopaedic procedures or other procedures where a thicker barrier* may be required
- Can be worn as a stand alone, when a single layer of protection is preferred or as the top or bottom gloves when double-gloving
- Synthetic, not made with natural rubber latex

Catalog No.	Size	Length	Thickness ¹			Material	Color	Cuff Type	Qty / bx	Qty / cs
			Size	Length	Finger					
2D73ET60	6	292mm	0.30mm	0.21mm	0.19mm	Synthetic polyisoprene (PI)	 Cream	Beaded/ Rolled	40	160
2D73ET65	6.5									
2D73ET70	7									
2D73ET75	7.5									
2D73ET80	8									
2D73ET85	8.5									
2D73ET80	9									



See Appendix page 29 for complete testing standards



*As compared to Protexis™ PI surgical gloves

1. Thickness tested in accordance with ASTM D 3577




CardinalHealth
Essential to care™



This product is not made
with natural rubber latex

Protexis™ Neoprene Surgical Gloves

- Thinner and softer for enhanced tactile response*
- Synthetic, not made with natural rubber latex

Catalog no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D73DP55	5.5	280mm	0.17mm	0.14mm	0.14mm	Synthetic neoprene with nitrile polymer coating	 Light brown	Beaded/ Rolled	50	200
2D73DP60	6									
2D73DP65	6.5									
2D73DP70	7	296mm								
2D73DP75	7.5									
2D73DP80	8									
2D73DP85	8.5									
2D73DP90	9									



See Appendix page 29 for complete testing standards



See Appendix page 28 for chemotherapy agent permeation results



*They are designed to be softer and thinner than previous Cardinal Health™ Neoprene Surgical Gloves, while maintaining physical property performance.

1. Thickness tested in accordance with ASTM D 3577

Protexis™ Surgical Gloves



CardinalHealth
Essential to care™



This product is not made with natural rubber latex

Protexis™ Neoprene Essential Surgical Gloves

- Manufactured without traditional chemical accelerators that have been known to cause skin sensitivities²
- Smooth finish for tactile sensitivity
- Nitrile coating for strength, protection and easy donnability
- Synthetic, not made with natural rubber latex

Catalog no.	Size	Length	Thickness ¹			Material	Color	Cuff type	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff					
2D73DS55	5.5	281mm	0.17mm	≥ 0.14mm	≥ 0.14 mm	Synthetic neoprene with nitrile polymer coating	 Light brown	Beaded/ Rolled	50	200
2D73DS60	6									
2D73DS65	6.5									
2D73DS70	7	295mm								
2D73DS75	7.5									
2D73DS80	8									
2D73DS85	8.5									
2D73DS90	9									



See Appendix page 29 for complete testing standards



Not made with chemical accelerators



See Appendix page 28 for chemotherapy agent permeation results



1. Thickness tested in accordance with ASTM D 3577

2. The chemical accelerators commonly used during traditional manufacturing of surgical gloves are DPG, MBTs, Thiurams, and Carbamates.




CardinalHealth
Essential to care™



NATURAL RUBBER LATEX¹

Protexis™ Latex Surgical Gloves

- Versatile glove to be used in a wide variety of surgical environments
- Brown tint blends with the wearer's skin for protection that is unobtrusive and less noticeable²

Catalog no.	Size	Length	Thickness ³			Material	Color	Cuff type	Protein content ⁴	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff						
2D72NS55X	5.5	282mm	0.25mm	0.20mm	0.19mm	Natural rubber latex with nitrile polymer coating	 Light Brown	Beaded/ Rolled	Less than 50 micrograms/ gram	50	200
2D72NS60X	6										
2D72NS65X	6.5										
2D72NS70X	7	295mm									
2D72NS75X	7.5										
2D72NS80X	8										
2D72NS85X	8.5										
2D72NS90X	9										



See Appendix page 29 for complete testing standards



See Appendix page 28 for chemotherapy agent permeation results



1. Caution: This product contains natural rubber latex which may cause allergic reactions, including anaphylactic response

2. As compared to other colored gloves in the Protexis™ Latex Surgical Gloves portfolio

3. Thickness tested in accordance with ASTM D 3577

4. Protein content tested in accordance with EN455-3 using ASTM D 5712




NATURAL RUBBER LATEX¹

Protexis™ Latex with Neu-Thera™ Surgical Gloves

- Designed to protect in a broad range of cases
- Same great engineering as our Protexis™ Latex Surgical Gloves with our patented Neu-Thera™ Coating

⑥ What is Neu-Thera™ Coating?²

We place moisturizing Neu-Thera™ Coating on the inside of Protexis™ Latex with Neu-Thera™ Surgical Gloves to **help promote overall skin wellness and support hand health**. It provides relief of dry, flaky skin and minimizes skin shedding.

Catalog no.	Size	Length	Thickness ³			Material	Color	Cuff type	Protein content ⁴	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff						
2D73TP55	5.5	281mm	0.25mm	0.20mm	0.19mm	Natural rubber latex with nitrile polymer and Neu-Thera™ Emollient Coating	 Light brown	Beaded/ Rolled	Less than 50 micrograms/ gram	50	200
2D73TP60	6										
2D73TP65	6.5										
2D73TP70	7	293mm									
2D73TP75	7.5										
2D73TP80	8										
2D73TP85	8.5										
2D73TP90	9										



See Appendix page 29 for complete testing standards



See Appendix page 28 for chemotherapy agent permeation results



Protexis™ Surgical Gloves

1. Caution: This product contains natural rubber latex which may cause allergic reactions, including anaphylactic response

2. Data on file

3. Thickness tested in accordance with ASTM D 3577

4. Protein content tested in accordance with EN455-3 using ASTM D 5712



CardinalHealth

Essential to care™



NATURAL RUBBER LATEX¹

Protexis™ Latex Essential Surgical Gloves

- Versatile glove to be used in a wide variety of surgical environments
- Excellent choice as an outer glove when double-gloving
- Thinner 20% than Protexis™ Latex Surgical Gloves for enhanced tactile sensitivity

Catalog no.	Size	Length	Thickness ²			Material	Color	Cuff type	Protein content ³	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff						
2D72LE55	5.5	280mm	0.19mm	0.21mm	0.16mm	Natural rubber latex with nitrile polymer coating		Beaded/ Rolled	Less than 50 micrograms/ gram	50	200
2D72LE60	6										
2D72LE65	6.5										
2D72LE70	7	292mm									
2D72LE75	7.5										
2D72LE80	8										
2D72LE85	8.5										
2D72LE90	9										



See Appendix page 29 for complete testing standards



See Appendix page 28 for chemotherapy agent permeation results



Protexis™ Surgical Gloves

1. Caution: This product contains natural rubber latex which may cause allergic reactions, including anaphylactic response

2. Thickness tested in accordance with ASTM D 3577

3. Protein content tested in accordance with EN455-3 using ASTM D 5712



CardinalHealth

Essential to care™




NATURAL RUBBER LATEX¹

Protexis™ Latex Blue with Neu-Thera™ Surgical Gloves

- Good underglove solution when double-gloving
- Distinct blue color aids in alerting wearers to perforations in the outer glove
- Designed to be comfortable and reliable for any type of surgical procedure

⑥ What is Neu-Thera™ Coating?²

We place moisturizing Neu-Thera™ Coating on the inside of Protexis™ PI Latex Blue with Neu-Thera™ Surgical Gloves to **help promote overall skin wellness and support hand health**. It provides relief of dry, flaky skin and minimizes skin shedding.

Catalog no.	Size	Length	Thickness ³			Material	Color	Cuff type	Protein content ⁴	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff						
2D72LU55	5.5	282mm	0.19mm	0.14mm	0.14mm	Natural rubber latex with nitrile polymer and Neu-Thera™ Emollient Coating	 Blue	Beaded/ Rolled	Less than 50 micrograms/gram	50	200
2D72LU60	6										
2D72LU65	6.5										
2D72LU70	7	295mm									
2D72LU75	7.5										
2D72LU80	8										
2D72LU85	8.5										
2D72LU90	9										

① See Appendix page 29 for complete testing standards



Protexis™ Surgical Gloves

1. Caution: This product contains natural rubber latex which may cause allergic reactions, including anaphylactic response

2. Data on file

3. Thickness tested in accordance with ASTM D 3577

4. Protein content tested in accordance with EN455-3 using ASTM D 5712


CardinalHealth
Essential to care™



NATURAL RUBBER LATEX¹

Protexis™ Latex Micro Surgical Gloves

- Stretches and conforms to your hand contour, staying comfortably in place
- Good solution in a thin double-gloving system where fingertip sensation is essential
- 30 percent thinner for enhanced flexibility and tactile sensitivity²
- Heightened tactile response with a comfortable, smooth, anti-slip finish

Catalog no.	Size	Length	Thickness ³			Material	Color	Cuff type	Protein content ⁴	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff						
2D72NT55X	5.5	282mm	0.17mm	0.14mm	0.14mm	Natural rubber latex with nitrile polymer coating	 Light brown	Beaded/ Rolled	Less than 50 micrograms/ gram	50	200
2D72NT60X	6										
2D72NT65X	6.5										
2D72NT70X	7	295mm									
2D72NT75X	7.5										
2D72NT80X	8										
2D72NT85X	8.5										
2D72NT90X	9										



See Appendix page 29 for complete testing standards



See Appendix page 28 for chemotherapy agent permeation results



1. Caution: This product contains natural rubber latex which may cause allergic reactions, including anaphylactic response

2. As compared to Protexis™ Latex Surgical Gloves

3. Thickness tested in accordance with ASTM D 3577

4. Protein content tested in accordance with EN455-3 using ASTM D 5712

Protexis™ Surgical Gloves



CardinalHealth


Essential to care™



NATURAL RUBBER LATEX¹

Protexis™ Latex Ortho Surgical Gloves

- Thickest glove in the latex portfolio
- Smooth finish for tactile sensitivity
- Water-based hydrogel coating for easy donning
- Durable for broaching and tactile for pinning
- Rich brown color reduces glare from surgical lighting

Catalog no.	Size	Length	Thickness ²			Material	Color	Cuff type	Protein content ³	Qty/ bx	Qty/ cs
			Finger	Palm	Cuff						
2D72LT60	6	282mm	0.34mm	0.24mm	0.21mm	Natural rubber latex with water-based hydrogel polymer coating	 Brown	Beaded/ Rolled	Less than 50 micrograms/ gram	40	160
2D72LT65	6.5										
2D72LT70	7										
2D72LT75	7.5										
2D72LT80	8										
2D72LT85	8.5										
2D72LT90	9										



See Appendix page 29 for complete testing standards



Protexis™ Surgical Gloves

1. Caution: This product contains natural rubber latex which may cause allergic reactions, including anaphylactic response

2. Thickness tested in accordance with ASTM D 3577

3. Protein content tested in accordance with EN455-3 using ASTM D 5712



CardinalHealth
Essential to care™

Environmental stewardship

Protexis™ Surgical Gloves Perfect Fit is a commitment to practicing sustainable manufacturing and packaging practices and reducing our environmental impact in the global supply chain. Our culture of safety and best practice is not lost on the global ecosystem in which we operate.



Protexis™ Surgical Gloves

SUSTAINABLE MANUFACTURING & PACKAGING PRACTICES¹

95% of surgical glove inner wallets are made from recyclable and renewable paper sources

Product shipping boxes utilize at least **80%** recycled corrugated cardboard

Up to **40%** of water used in Cardinal Health surgical gloves manufacturing is recycled and reused on future production runs

All Cardinal Health™ Protexis™ Surgical Gloves are made **DEHP and PVC-free**

Cardinal Health was the first surgical gloves manufacturer to launch **sustainable half-fold packaging**



Cardinal Health™ Protexis™ Surgical Gloves

Testing standards and technical data

Non-latex polyisoprene



This product is not made
with natural rubber latex



PROPERTIES (BEFORE AGING)

	Protexis™ PI	Protexis™ PI with Neu-Thera™	Protexis™ PI Blue with Neu-Thera™	Protexis™ PI Micro	Protexis™ PI Orthopaedic	Protexis™ PI Ortho
Force at break	19N ¹	19N ¹	19N ¹	17N ¹	29N ¹	26N ¹
Elongation at break (elasticity) (min)	650% ²	650% ²	650% ²	650% ²	650% ²	650% ²
Puncture resistance (cuff) (min)	5N ³	5N ³	5N ³	5N ³	5N ³	5N ³
Freedom from holes	0.65 AQL ⁴	0.65 AQL ⁴	0.65 AQL ⁴	0.65 AQL ⁴	0.65 AQL ⁴	0.65 AQL ⁴
Sterilization	Gamma radiation	Gamma radiation	Gamma radiation	Gamma radiation	Gamma radiation	Gamma radiation
Accelerant	Zinc diethyldithiocarbamate (ZDEC), Zinc mercaptobenzothiazole (ZMBT), Diphenylguanidine (DPG)					

Non-latex neoprene



This product is not made
with natural rubber latex



PROPERTIES (BEFORE AGING)

	Protexis™ Neoprene	Protexis™ Neoprene Essential
Force at break	11N ¹	13N ¹
Elongation at break (elasticity) (min)	650% ²	650% ²
Puncture resistance (cuff) (min)	5N ³	—
Freedom from holes	0.65 AQL ⁴	0.65 AQL ⁴
Sterilization	Gamma radiation	Gamma radiation
Accelerant	ZDBC (Zinc Dibutyldithiocarbamate)	Not made with traditional chemical accelerators

1. Force at break tested in accordance with EN455-2
2. Elongation at break tested in accordance with ASTM D 3577
3. Puncture resistance (cuff) tested in accordance with AS/NZS 4179
4. Freedom from holes test in accordance with EN455-1

Latex



NATURAL RUBBER LATEX⁵



PROPERTIES (BEFORE AGING)

	Protexis™ Latex	Protexis™ Latex with NeuThera™	Protexis™ Latex Essential	Protexis™ Latex Blue with NeuThera™	Protexis™ Latex Micro	Protexis™ Latex Ortho
Force at break	20N ¹	19N ¹	17N ¹	18N ¹	17N ¹	34N ¹
Elongation at break (elasticity) (min)	750% ²	650% ²	750% ²	750% ²	750% ²	750% ²
Puncture resistance (cuff) (min)	5N ³	5N ³	9N ³	5N ³	5N ³	5N ³
Freedom from holes	0.65 AQL ⁴	0.65 AQL ⁴	0.65 AQL ⁴	0.65 AQL ⁴	0.65 AQL ⁴	0.65 AQL ⁴
Sterilization	Gamma radiation	Gamma radiation	Gamma radiation	Gamma radiation	Gamma radiation	Gamma radiation
Accelerant	ZDBC (Zinc Dibutyldithiocarbamate)					

1. Force at break tested in accordance with EN455-2

2. Elongation at break tested in accordance with ASTM D 3577

3. Puncture resistance (cuff) tested in accordance with AS/NZS 4179

4. Freedom from holes test in accordance with EN455-1

5. Caution: This product contains natural rubber latex which may cause allergic reactions, including anaphylactic response





Chemotherapy agent permeation testing

Agent	Minimum breakthrough detection time in minutes (0.01 µg/cm ² /minute) ¹							
	Protexis™ PI	Protexis™ PI Micro	Protexis™ Neoprene	Protexis™ Neoprene Essential	Protexis™ Latex	Protexis™ Latex Micro	Protexis™ Latex with Neu-Thera™	Protexis™ Latex Essential
Carmustine (3.3 mg/mL)	15.26 (Do Not Use)	15.5 (Do Not Use)	31.1 (Do Not Use)	60.1	15.2 (Do Not Use)	12.5 (Do Not Use)	15.0 (Do Not Use)	23.0 (Do Not Use)
Cisplatin (1.0 mg/mL)	>240	>240	>240	>240	>240	>240	>240	>240
Cyclophosphamide (20 mg/mL)	>240	>240	>240	>240	>240	>240	>240	>240
Doxorubicin HCL (2.0 mg/mL)	>240	>240	>240	>240	>240	>240	>240	>240
Etoposide (20 mg/mL)	>240	>240	>240	>240	>240	>240	>240	>240
Fluorouracil (50 mg/mL)	>240	>240	>240	>240	>240	>240	>240	>240
Ifosfamide (50 mg/mL)	>240	>240	>240	>240	>240	>240	>240	>240
Methotrexate (25 mg/mL)	>240	>240	>240	>240	>240	>240	>240	>240
Mitomycin C (0.5 mg/mL)	>240	>240	>240	>240	>240	>240	>240	>240
Mitoxantrone (2.0 mg/mL)	>240	>240	>240	>240	>240	>240	>240	>240
Paclitaxel (6.0 mg/mL)	>240	>240	>240	>240	>240	>240	>240	>240
Thiotepa (10 mg/mL)	16.04 (Do Not Use)	15.7 (Do Not Use)	76.0	110.5	16.3 (Do Not Use)	22.9 (Do Not Use)	20.5 (Do Not Use)	>240
Vincristine Sulfate (1.0 mg/mL)	>240	>240	>240	>240	>240	>240	>240	>240



When chemotherapy drugs are present, glove selection should be based on the specific type(s) of chemicals used. Users must refer to the instructions for use supplied with the box, for the chemicals being used to determine an adequate level of protection.

1. These gloves have been tested for resistance to permeation of various chemotherapy drugs per ASTM D 6978, "Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs."

Testing standards

GLOBAL QUALITY STANDARDS TESTED FOR AND ADHERED TO (results on file)

ASTM D3577, EN 455-2, ISO 10282	Physical dimension (length, width, palm)
ASTM D3577, EN 455-2, ISO 10282	Physical properties (tensile strength)
ASTM D624	Tear strength testing (T-tear, V-tear)
AS/NZA 4179	Puncture resistant testing
ASTM D5151, EN 455-1, ISO 10282	Freedom from holes (water-tightness)
ASTM D6124, EN ISO 21171	Powder residue for powder-free gloves
ASTM D6124	Powder amount for powdered gloves
ASTM D5712	Aqueous extractable protein content
ASTM D6499	Antigenic protein content
EN 455-3	Leachable protein level, modified Lowry method
ASTM D7102, EN 455-3	Endotoxin
ASTM D7160	Storage stability, accelerated aging
ASTM D7161	Storage stability, real-time aging
ASTM F739, EN 16523	Lab chemical permeation
ASTM D6978, EN 16523	Chemotherapy drug permeation (results on page 28)
ASTM F1671	Bacteriophage penetration
ISO 10993-10	Sensitivity testing & primary skin irritation
ISO 10993-5	Cytotoxicity testing
ISO 10993-11	Systemic Toxicity Testing (Acute)
ASTM D3577	Sterility test

PPE REQUIRED TESTING (results on file)

EN 420:2003 + A1:2009	General requirements, size, dexterity, pH and extractable protein
EN 388:2016	Protective gloves against mechanical risks
EN ISO 374-1:2016	Performance requirements for chemicals risk
EN 374-2:2014	Resistance to penetration against dangerous chemicals and micro-organisms
EN 16523-1:2015	Materials resistance to permeation by chemicals
EN 374-4:2013	Resistance to degradation by chemicals
EN ISO 374-5:2016	Performance requirements for micro-organisms risk

INTERNAL PROCEDURES (data on file)







- Bone cement permeation
- Low-hydration conductivity
- Residue accelerator test



Non-latex polyisoprene



This product is not made with natural rubber latex






	Catalog no.	Size	Length	Thickness*			Material <i>Coating in red</i>	Color	Cuff type	Chemical accelerant	Qty/ bx	Qty/ cs	
				Finger	Palm	Cuff							
	Protexis™ PI	2D72PT55X	5.5	287mm	0.23mm	0.17mm	0.17mm	Synthetic polyisoprene (PI)	 Cream	Beaded/ Rolled	1. 1, 3-Diphenylguanidine (DPG) 2. Zinc Diethyldithiocarbamate (ZDEC) 3. Zinc-2-mercaptobenzothiazole (ZMBT)	50	200
		2D72PT60X	6										
		2D72PT65X	6.5										
		2D72PT70X	7	300mm									
		2D72PT75X	7.5										
		2D72PT80X	8										
		2D72PT85X	8.5										
		2D72PT90X	9										
	Protexis™ PI with Neu-Thera™	2D73TE55	5.5	287mm	0.23mm	0.25mm	0.17mm	Synthetic polyisoprene (PI) with Neu-Thera™ Emollient Coating	 Cream	Beaded/ Rolled	1. 1, 3-Diphenylguanidine (DPG) 2. Zinc Diethyldithiocarbamate (ZDEC) 3. Zinc-2-mercaptobenzothiazole (ZMBT)	50	200
		2D73TE60	6										
		2D73TE65	6.5										
		2D73TE70	7	300mm									
		2D73TE75	7.5										
		2D73TE80	8										
		2D73TE85	8.5										
		2D73TE90	9										
	Protexis™ PI Blue with Neu-Thera™	2D73EB55	5.5	287mm	0.20mm	0.14mm	0.14mm	Synthetic polyisoprene (PI) with Neu-Thera™ Emollient Coating	 Blue	Beaded/ Rolled	1. 1, 3-Diphenylguanidine (DPG) 2. Zinc Diethyldithiocarbamate (ZDEC) 3. Zinc-2-mercaptobenzothiazole (ZMBT)	50	200
		2D73EB60	6										
		2D73EB65	6.5										
		2D73EB70	7	300mm									
		2D73EB75	7.5										
		2D73EB80	8										
		2D73EB85	8.5										
		2D73EB90	9										

*Thickness tested in accordance with ASTM D 3577

Non-latex polyisoprene



This product is not made with natural rubber latex

	Catalog no.	Size	Length	Thickness*			Material Coating in red	Color	Cuff type	Chemical accelerant	Qty/ bx	Qty/ cs	
				Finger	Palm	Cuff							
	Protexis™ PI Micro	2D73PM55	5.5	287mm	0.20mm	0.14mm	0.14mm	Synthetic polyisoprene (PI)	 Cream	Beaded/ Rolled	1. 1, 3-Diphenylguanidine (DPG) 2. Zinc Diethyldithiocarbamate (ZDEC) 3. Zinc-2-mercaptobenzothiazole (ZMBT)	50	200
		2D73PM60	6										
		2D73PM65	6.5										
		2D73PM70	7	300mm									
		2D73PM75	7.5										
		2D73PM80	8										
		2D73PM85	8.5										
		2D73PM90	9										
	Protexis™ PI Orthopaedic	2D73HT60	6	291mm	0.34mm	0.26mm	0.21mm	Synthetic polyisoprene (PI) with water-based hydrogel polymer coating	 Brown	Beaded/ Rolled	1. 1, 3-Diphenylguanidine (DPG) 2. Zinc Diethyldithiocarbamate (ZDEC) 3. Zinc-2-mercaptobenzothiazole (ZMBT)	160	200
		2D73HT65	6.5										
		2D73HT70	7	305mm									
		2D73HT75	7.5										
		2D73HT80	8										
		2D73HT85	8.5										
		2D73HT90	9										
			Protexis™ PI Ortho	2D73ET60									
2D73ET65	6.5												
2D73ET70	7			305mm									
2D73ET75	7.5												
2D73ET80	8												
2D73ET85	8.5												
2D73ET90	9												



*Thickness tested in accordance with ASTM D 3577

Non-latex neoprene



This product is not made with natural rubber latex









	Catalog no.	Size	Length	Thickness*			Material <i>Coating in red</i>	Color	Cuff type	Chemical accelerant	Qty/ bx	Qty/ cs
				Finger	Palm	Cuff						
Protexis™ Neoprene	2D73DP55	5.5	280mm	0.17mm	0.14mm	0.14mm	Synthetic neoprene with nitrile polymer	 Light Brown	Beaded/ Rolled	Zinc Dibutyldithiocarbamate (ZDBC)	50	200
	2D73DP60	6										
	2D73DP65	6.5										
	2D73DP70	7	296mm									
	2D73DP75	7.5										
	2D73DP80	8										
	2D73DP85	8.5										
	2D73DP90	9										
Protexis™ Neoprene Essential	2D73DS55	5.5	281mm	0.17mm	≥ 0.14mm	≥ 0.14 mm	Synthetic neoprene with nitrile polymer	 Light Brown	Beaded/ Rolled	Not made with traditional chemical accelerators	50	200
	2D73DS60	6										
	2D73DS65	6.5										
	2D73DS70	7	295mm									
	2D73DS75	7.5										
	2D73DS80	8										
	2D73DS85	8.5										
	2D73DS90	9										

*Thickness tested in accordance with ASTM D 3577

Latex









	Catalog no.	Size	Length	Thickness*			Material Coating in red	Color	Cuff type	Protein content	Chemical accelerant	Qty/ bx	Qty/ cs	
				Finger	Palm	Cuff								
	Protexis™ Latex	2D72NS55X	5.5	282mm	0.25mm	0.20mm	0.19mm	Natural rubber latex with nitrile polymer	 Light Brown	Beaded/ Rolled	Less than 50 micrograms/ gram	Zinc Dibutylthio- carbamate (ZDBC)	50	200
		2D72NS60X	6											
		2D72NS65X	6.5											
		2D72NS70X	7	295mm										
		2D72NS75X	7.5											
		2D72NS80X	8											
		2D72NS85X	8.5											
		2D72NS90X	9											
	Protexis™ Latex with Neu-Thera™	2D73TP55	5.5	281mm	0.25mm	0.20mm	0.19mm	Natural rubber latex with nitrile polymer and Neu-Thera™ Emollient Coating	 Light Brown	Beaded/ Rolled	Less than 50 micrograms/ gram	Zinc Dibutylthio- carbamate (ZDBC)	50	200
		2D73TP60	6											
		2D73TP65	6.5											
		2D73TP70	7	293mm										
		2D73TP75	7.5											
		2D73TP80	8											
		2D73TP85	8.5											
		2D73TP90	9											
	Protexis™ Latex Essential	2D72LE55	5.5	280mm	0.19mm	0.21mm	0.16mm	Natural rubber latex with nitrile polymer	 Cream	Beaded/ Rolled	Less than 50 micrograms/ gram	Zinc Dibutylthio- carbamate (ZDBC)	50	200
		2D72LE60	6											
		2D72LE65	6.5											
		2D72LE70	7	292mm										
		2D72LE75	7.5											
		2D72LE80	8											
		2D72LE85	8.5											
		2D72LE90	9											

*Thickness tested in accordance with ASTM D 3577

Latex



	Catalog no.	Size	Length	Thickness*			Material Coating in red	Color	Cuff type	Protein content	Chemical accelerant	Qty/ bx	Qty/ cs	
				Finger	Palm	Cuff								
	Protexis™ Latex Blue with Neu-Thera™	2D72LU55	5.5	282mm	0.19mm	0.14mm	0.14mm	Natural rubber latex with nitrile polymer and Neu-Thera™ Emollient Coating	 Blue	Beaded/ Rolled	Less than 50 micrograms/ gram	Zinc Dibutyldithio-carbamate (ZDBC)	50	200
		2D72LU60	6											
		2D72LU65	6.5											
		2D72LU70	7	295mm										
		2D72LU75	7.5											
		2D72LU80	8											
		2D72LU85	8.5											
		2D72LU90	9											
	Protexis™ Latex Micro	2D72NT55X	5.5	282mm	0.17mm	0.14mm	0.14mm	Natural rubber latex with nitrile polymer	 Light Brown	Beaded/ Rolled	Less than 50 micrograms/ gram	Zinc Dibutyldithio-carbamate (ZDBC)	50	200
		2D72NT60X	6											
		2D72NT65X	6.5											
		2D72NT70X	7	295mm										
		2D72NT75X	7.5											
		2D72NT80X	8											
		2D72NT85X	8.5											
		2D72NT90X	9											
	Protexis™ Latex Ortho	2D72LT60	6	282mm	0.34mm	0.24mm	0.21mm	Natural rubber latex with water-based hydrogel polymer coating	 Brown	Beaded/ Rolled	Less than 50 micrograms/ gram	Zinc Dibutyldithio-carbamate (ZDBC)	40	160
		2D72LT65	6.5											
		2D72LT70	7	295mm										
		2D72LT75	7.5											
		2D72LT80	8											
		2D72LT85	8.5											
		2D72LT90	9											

*Thickness tested in accordance with ASTM D 3577



QUALITY

Protexis™ Surgical Gloves

Perfect Fit



CLINICAL BEST PRACTICE



TRUST

For healthcare professionals only

Important information: Prior to use, refer to the instructions on the dispenser box supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cardinal Health reserves the right to change product specifications without prior notification.

Please contact your Cardinal Health representative for additional product availability information.

Not for distribution in the United States.

© 2018 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, ESSENTIAL TO CARE, PROTEXIS and NEU-THERA are trademarks of Cardinal Health and may be registered in the US and/or in other countries. Lit. No. 2MP18-878870 (01/2019)