

620100 Surgical Hood Tuck

Surgical Hood

Product details

Size :	One size
Descriptive feature :	Standard, Tie Band
Colour :	Green
Sterile :	Non-sterile
Brand :	BARRIER

Images



Delivered items

620100-40

Sales released in: Algeria, Austria, Belgium, Bulgaria, Cyprus, Czechia, Denmark, Estonia, Finland, France, French Polynesia, Germany, Greece, Guadeloupe, Hong Kong, Hungary, Italy, Kuwait, Latvia, Lithuania, Luxembourg, Martinique, Mauritius, Moldova (the Republic of), Monaco, Netherlands, New Caledonia, Norway, Oman, Poland, Portugal, Qatar, Romania, Singapore, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey, United Arab Emirates, United Kingdom of Great Britain and Northern Ireland

Country of origin: Morocco

Shelf life: 5 years

Sterilization method: Non-sterile

Production Responsibility: Lamatem, Ard Lakhloualet. Zone Industrielle Berrechid, Morocco 26300

Packing information: First packaging layer is a cardboard dispenser box. Second layer is a corrugated board transport box.

Is suitable for Tray: No

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2022 Mölnlycke Health Care AB. All rights reserved.



Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Piece	1	7323190257010			
Dispenser box	120	7323190221219			
Transport box	960	7323190221202	265x593x439 mm	69.0 dm3	6.3 / 4.8 kg
Pallet	23040	7323190221196	800x1200x1906 mm		176.7 / 115.2 kg

Material

Animal tissues :	No
Human blood derivatives :	No
Natural rubber latex :	No
Medicinal substances :	No
Phthalates :	No
Polyvinyl chloride :	No

Product Composition Headwear

Product Component	Composition
Main material	Viscose nonwoven, 25 gsm
Sewing thread	Polyester

Technical

Dimension

Dimension reference	Dimension text	Dimension value
A	Length straps	62 cm
B	Length side piece	28 cm
C	Face length	17 cm
D	Face width	9,5 cm
E	Width tie-band	4 cm
F	Length top piece	25 cm
G	Width top piece	15 cm
H	Height forehead	4 cm

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2022 Mölnlycke Health Care AB. All rights reserved.



Disposal instructions

Non-hazardous waste used BARRIER products and sterility barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. BARRIER products do not contain any hazardous substances that can leach out if the products are land filled. Transport boxes are designed to fit existing recovery systems. The new BARRIER packaging system complies with the Packaging Waste Directive of the European Union.

Storage instructions

Mölnlycke Health Care recommends that BARRIER products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterility barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the package is opened and the product is used.

Classifications

Regulation type(s)	MDR Class I ns	Locally Regulated	Locally Regulated	Locally Regulated
Intended Purpose :	The intended use of headwear is infection control by reducing the risk of transferring bacteria, hair or particles into the surgical wound	The intended use of headwear is infection control by reducing the risk of transferring bacteria, hair or particles into the surgical wound	The intended use of headwear is infection control by reducing the risk of transferring bacteria, hair or particles into the surgical wound	The intended use of headwear is infection control by reducing the risk of transferring bacteria, hair or particles into the surgical wound
Notified body medical devices/PPE :	MPA	Local Authority Kuwait	Local Authority MEA	Local Authority MEA
510(k) clearance number :			NA	
Regulatory Released :	Austria, Belgium, Bulgaria, Cyprus, Czechia, Denmark, Estonia, Finland, France, Germany, Greece, Guadeloupe, Hungary, Italy, Latvia, Lithuania, Luxembourg, Martinique, Monaco, Netherlands, New Caledonia, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland, United Kingdom of Great Britain and Northern Ireland	Kuwait	Turkey	Algeria, Mauritius, Oman, Qatar, Turkey

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.

Classifications

Regulation type(s)	Locally Regulated	Locally Regulated	Locally Regulated
Intended Purpose :	The intended use of headwear is infection control by reducing the risk of transferring bacteria, hair or particles into the surgical wound	The intended use of headwear is infection control by reducing the risk of transferring bacteria, hair or particles into the surgical wound	The intended use of headwear is infection control by reducing the risk of transferring bacteria, hair or particles into the surgical wound
Notified body medical devices/PPE :	Local Authority United Arab Emirates	Local Authority APAC	Local Authority CIS
510(k) clearance number :			
Regulatory Released :	United Arab Emirates	French Polynesia, Hong Kong, Singapore	Moldova (the Republic of)

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 10993-5, EN ISO 15223-1, ISO 15223-2, EN ISO 10993-10, ISO 14001

GMDN Code (Global Medical Device Nomenclature)

32297 Surgical cap, single-use

UNSPSC

42131611 Surgeon caps or hoods

Commodity Code

6505009090 Caps and hoods

CE Responsibility

Mölnlycke Health Care AB, Gamlestadsvägen 3c, Box 13080, SE-402 52 Goteborg Sweden

Basic UDID

733243000000000207JY

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.
©2022 Mölnlycke Health Care AB. All rights reserved.

