BD Discardit[™] II syringe

The BD Discardit" II syringe is a single-use medical device intended for injection and aspiration of medical fluids, including both bodily fluids (blood, etc.) and medicines. The lubricant included in the barrel material allows for smooth advancement of the plunger, leading to reduced sliding force.

2-Piece Design

The robust retaining ring helps to reduce the incidence of accidental plunger rod pullout during aspiration while the clear and transparent barrel provides easy visualization of the syringe content and dosing accuracy.

- Transparent barrel for good visualization of syringe content
- Smooth plunger movement, reduced sliding force
- Leak-tight
- Accurate graduation suited to user's needs
- Retaining ring to prevent accidental plunger rod withdrawal
- Design of finger grip flanges ensures stability and comfort during injection
- Ergonomically designed plunger for single-handed operation
- Sterile, single-use
- Natural Rubber Latex is not part of the material formulation



BD Discardit" II Syringes									
Reference	Capacity (ml)	Tip	Scale (ml)	Units / Box					
300928	2	Luer slip concentric	0.1	100					
309050	5	Luer slip eccentric	0.2	100					
309110	10	Luer slip eccentric	0.5	100					
300296	20	Luer slip eccentric	1	80					

BD Discardit" II Syringes with detached BD Microlance" Needle											
Reference	Gauge	Capacity (ml)	Tip	Length (mm)	Length (inch)	Colour code	Scale (ml)	Units / Box			
300322	23 G	2	Luer slip concentric	25	1"	blue	0.1	100			
300320	22 G	2	Luer slip concentric	30	1 ¼"	black	0.1	100			
300318	21 G	2	Luer slip concentric	40	1 ½"	green	0.1	100			
300325	21 G	2	Luer slip concentric	16	5/8"	green	0.1	100			
300516	23 G	5	Luer slip eccentric	30	1 ¼"	blue	0.2	100			
300513	22 G	5	Luer slip eccentric	30	1 ¼"	black	0.2	100			
300512	21 G	5	Luer slip eccentric	40	1 ½"	green	0.2	100			
300626	22 G	10	Luer slip eccentric	30	1 ¼"	black	0.5	100			
301001	21 G	10	Luer slip eccentric	40	1 ½"	green	0.5	100			
300330	21 G	20	Luer slip eccentric	40	1 ½"	green	1	80			

